



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stand ring 70 mm ø

€9.12
plus VAT &
Shipping

Product Images



Description

The **ecoLab stand rings** made of powder-coated steel are robust and ideal for hanging separating funnels and other conical devices. Select the ring with the appropriate diameter for your needs. This allows you to separate phases without worry over an extended period of time.

- 12 mm stem diameter

Additional Information

| | |
|-----------------------|-------------------|
| No. | E-0132 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072121082 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Steel galvanized |
| Form | annular. closed |
| suitable for | Separating funnel |
| MAX span | 70 mm |
| Span MIN | 70 mm |
| for DM Tripod Rod MAX | 12 mm |
| for DM tripod rod MIN | 8 mm |
| for clamping to | Tripod rods |

