



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Usbeck tripod rod stainless steel 75 cm long

Special Price

€12.60 was

€16.40

**plus VAT &
Shipping**

Product Images



Description

Tripod rods made of 18/8 stainless steel
With M10 thread.
Rod diameter 12 mm.

Additional Information

No.	UB-0116
Manufacturer (Brand)	Usbeck
old neolab article no.	2-2296
EAN	4058072399191
Transport temperature	Room temperature
Material	Stainless steel
Thread	M 10
DM Rod	12 mm
Length rod / arm	750 mm
Type rod	Tripod rod

