




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoFroxx Peptide blocking solution for immunodetection

€169.60 As low
as
€150.20
plus VAT &
Shipping

Product Images

Description

Blockers based on chemically modified serum proteins; powerful alternative to BSA blocking solutions.

Western blot membranes must bind proteins and antibodies optimally so that even low protein concentrations can be detected effectively. However, due to the binding property of the membrane surface, insufficient blocking of the free binding sites after protein transfer results in non-specific background signals and erroneous results. For clean blots and reproducible, reliable results, suitable blocking is therefore mandatory. Depending on the goal of the application (routine; quality control; publication; diagnostics) and the type of sample, different types of blockers are available. For many standard applications, the inexpensive milk powder or albumin (BSA) is the blocker of choice. However, for high background and non-specific bands, a switch to more potent alternatives is recommended.

Additional Information

| | |
|-----------------------|------------------|
| No. | 1095ML500 |
| Manufacturer (Brand) | neoFroxx |
| EAN | 4066306026980 |
| HSNumber | 38220000 |
| autoclavable at 121°C | No |
| Pack | 500 ml |
| Storage | Room temperature |
| sterile | No |

