



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoFroxx WesternFroxx all-in-one protein marker (15 - 200 kDa) for immunodetection

**€243.30** As low  
as  
**€215.50**  
plus VAT &  
Shipping

### Product Images

## Description

With the **WesternFroxx all-in-one protein marker (15 – 200 kDa)** we offer you a protein size standard that is not only pre-stained but also detectable via antibody detection analogous to the target protein. The result: clear band assignment; no more inaccuracies due to manual labeling.

A must for the perfect Western blot! The advantages at a glance:

- Prestained: 4 colored bands for visual control on gel and membrane
- Immunodetectable: 10 IgG binding proteins allow visualization of the all-in-one marker with the respective detection method used (IgG binding sites are used to bind the primary/secondary antibodies to the marker proteins during antibody incubation)
- Compatible with ECL-, fluorescence- and chromophore-based detection systems
- No additional antibodies/reagents/work steps
- Suitable for stripping & Reprobing
- Ready-to-use (application volume 1 – 2.5 µl )

Our Tip:

Combine the all-in-one protein marker with our WesternFroxx kits (article no. 5555 and 6666) for nicer and faster blots

## Additional Information

No.	7777YL250
neoGreen	No
Manufacturer (Brand)	neoFroxx
EAN	4066306031717
HSNumber	38220000
Transport temperature	Wet ice
autoclavable at 121°C	No
Pack	250 µl
Storage	Room temperature
sterile	No

