



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoFroxx WesternFroxx all-in-one protein marker (15 - 200 kDa) for immunodetection

€243.30 As low
as
€215.50
plus VAT &
Shipping

Product Images

Description

With the **WesternFroxx all-in-one protein marker (15 - 200 kDa)** we offer you a protein size standard that is not only pre-stained but also detectable via antibody detection analogous to the target protein. The result: clear band assignment; no more inaccuracies due to manual labeling.

A must for the perfect Western blot! The advantages at a glance:

- Prestained: 4 colored bands for visual control on gel and membrane
- Immunodetectable: 10 IgG binding proteins allow visualization of the all-in-one marker with the respective detection method used (IgG binding sites are used to bind the primary/secondary antibodies to the marker proteins during antibody incubation)
- Compatible with ECL-, fluorescence- and chromophore-based detection systems
- No additional antibodies/reagents/work steps
- Suitable for stripping & Reprobing
- Ready-to-use (application volume 1 - 2.5 µl)

Our Tip:

Combine the all-in-one protein marker with our WesternFroxx kits (article no. 5555 and 6666) for nicer and faster blots

Additional Information

No.	7777YL250
Manufacturer (Brand)	neoFroxx
EAN	4066306031717
Transport temperature	Wet ice
HSNumber	38220000
autoclavable at 121°C	No
Pack	250 µl
Storage	Room temperature
sterile	No

