




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Biological Industries EZ-PCR Mycoplasma Test Kit, 20 Tst

€277.60
plus VAT &
Shipping

Product Images

Description

Detection of Mycoplasma by PCR is nowadays the preferred method to detect the bacterium quickly, specifically and reliably in cell cultures and media. The EZ-PCR Mycoplasma Test Kit from Biological Industries meets all the requirements for regular routine testing for Mycoplasma: It is user-friendly (only 3 components: Ready-to-use PCR reaction mix with Taq polymerase, buffer solution for sample collection, and positive control), sensitive, specific, and provides a reliable result in less than 5 hours. The highly specific, 16S rRNA-based primers detect over 20 different mycoplasma species characteristic of cell culture contamination². The specific buffer composition and optimized sample preparation protocol also prevent PCR inhibition by media components, eliminating the need for additional internal control.

Additional Information

No.	20-700-20
Manufacturer (Brand)	Biological Industries
EAN	4058072983659
HSNumber	38220000
Transport temperature	Dry ice
autoclavable at 121°C	No
Pack	20 Tst
Storage	Room temperature
sterile	No

