



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Biological Industries EZ-PCR Mycoplasma Test Kit, 20 Tst

€1,298.00
plus VAT &
Shipping

Product Images

Description

The EZ-PCR Mycoplasma Kit is a user-friendly and convenient solution for routine mycoplasma detection in cell cultures. The test is based on amplification of a mycoplasma-specific 16S rRNA gene region and provides reliable and reproducible results. NEW: To indicate possible PCR inhibition by sample components, the kit now also includes an internal control DNA. This makes false-negative results a thing of the past!

- Ready-to-use: Minimization of pipetting errors
- Fast and easy to use: results after just a few hours
- Sensitive (10CFU/ml)
- Positive control (270 bp)
- Internal control (357 bp): to detect an intact PCR reaction

Additional Information

No.	20-700-20BX5
Manufacturer (Brand)	Biological Industries
Transport temperature	Dry ice
HSNumber	38220000
autoclavable at 121°C	No
Pack	5 x 20 Tst
Storage	Room temperature
sterile	No

