



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Acrylamide/Bisacrylamide Solution 40% (29:1), 100 ml

€22.00
plus VAT &
Shipping

Product Images

Description

Hi-Media Acrylamide/Bisacrylamide Solution 40% (29:1), 100 ml

Additional Information

| | |
|-----------------------|------------------|
| No. | ML084-100ML |
| Manufacturer (Brand) | Hi-Media |
| HSNumber | 38210000 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Pack | 100 ml |
| Storage | Cooling |
| sterile | No |

