



## HiMedia Adenine, 100 g

### Product Images

---



## Description

Hi-Media Adenine, 100 g

## Additional Information

|                      |                |
|----------------------|----------------|
| No.                  | TC082-100G     |
| ClassPG              | 6.1/III        |
| Manufacturer (Brand) | Hi-Media       |
| CAS number           | 73-24-5        |
| HSNumber             | 29335995       |
| Pack                 | 100 g          |
| Storage              | Raumtemperatur |
| UN number            | UN2811         |
| sterile              | No             |

