



## HiMedia 1,4-Dibromobutane, 250 ml

### Product Images

---



## Description

Hi-Media 1,4-Dibromobutane, 250 ml

## Additional Information

|                      |                |
|----------------------|----------------|
| No.                  | RM9431-250ML   |
| ClassPG              | 6.1            |
| Manufacturer (Brand) | Hi-Media       |
| HSNumber             | 29033919       |
| Pack                 | 250 ml         |
| Storage              | Raumtemperatur |
| UN number            | UN2810         |

