



## HiMedia DL-Asparagine, monohydrate, 25 g

### Product Images

---



## Description

Hi-Media DL-Asparagine, monohydrate, 25 g

## Additional Information

|                      |                |
|----------------------|----------------|
| No.                  | RM1475-25G     |
| Manufacturer (Brand) | Hi-Media       |
| CAS number           | 3130-87-8      |
| HSNumber             | 29241900       |
| Pack                 | 25 g           |
| Storage              | Raumtemperatur |
| sterile              | No             |

