



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia DL-Lactic acid Sodium salt solution, 60% Sterile filtered, 100 ml

€70.00
plus VAT &
Shipping

Product Images

Description

Hi-Media DL-Lactic acid Sodium salt solution, 60% Sterile filtered, 100 ml

Additional Information

No.	TCL138-100ML
Manufacturer (Brand)	Hi-Media
Transport temperature	Room temperature
autoclavable at 121°C	No
Pack	100 ml
Storage	Room temperature
sterile	Yes

