



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Fixative (BFA), for rapid fixing of haematological samples, 500 ml

€61.30
plus VAT &
Shipping

Product Images



Description

Hi-Media Fixative (BFA), for rapid fixing of haematological samples, 500 ml

Additional Information

| | |
|-----------------------|------------------|
| No. | S119-500ML |
| Manufacturer (Brand) | Hi-Media |
| HSNumber | 38220000 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Pack | 500 ml |
| Storage | Room temperature |
| sterile | No |

