



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Fixative (Buffered Formalin fixative) for fixing cytological or histolo

€56.90
plus VAT &
Shipping

Product Images



Description

Hi-Media Fixative (Buffered Formalin fixative) for fixing cytological or histological samples, 500 ml

Additional Information

No.	S109-500ML
Manufacturer (Brand)	Hi-Media
HSNumber	38220000
Transport temperature	Room temperature
autoclavable at 121°C	No
Pack	500 ml
Storage	Room temperature
sterile	No

