



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Geneticin Solution (G418 Sol.) w/ 50 mg/ml G418 in tissue culture grade

€230.30
plus VAT &
Shipping

Product Images

Description

Hi-Media Geneticin Solution (G418 Sol.) w/ 50 mg/ml G418 in tissue culture grade water, sterile, 20 ml

Additional Information

No.	A016-20ML
Manufacturer (Brand)	Hi-Media
Transport temperature	Dry ice
autoclavable at 121°C	No
Pack	20 ml
Storage	Room temperature
sterile	No

