



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia HiCrome™ Coliform Agar w/ SLS, 500 g

€358.80
plus VAT &
Shipping

Product Images

Description

Nutrient media for drinking water testing should cover as wide a range of microorganisms as possible, i.e. provide an optimum range of nutrients, in order to be able to cultivate as many germs as possible. According to the German Drinking Water Ordinance, a limit of 100 cfu/ml (colony-forming units per milliliter) applies to piped drinking water. Coliform germs (indicator of fecal contamination) have a special status; here, not a single colony may be detected per 100 ml.

Additional Information

No.	MI300-500G
Manufacturer (Brand)	Hi-Media
HSNumber	38210000
Transport temperature	Room temperature
autoclavable at 121°C	No
Pack	500 g
Storage	Cooling
sterile	No

