



HiMedia HiVeg™ Extract, 500 g

Product Images



Description

The starting material for all HiVeg™ products are selected plant proteins, which are digested to peptides and amino acids using plant enzymes. For each animal product, many years of work have gone into developing a plant-based substitute which is virtually indistinguishable from the original in terms of both composition and growth characteristics.

Additional Information

No.	RM002V-500G
Manufacturer (Brand)	Hi-Media
HSNumber	35040090
Pack	500 g
Storage	Raumtemperatur
sterile	No

