



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## HiMedia Pseudomonas Agar For Detection of Pyocyanin, 500 g

**Price on  
request  
plus VAT &  
Shipping**

### Product Images

## Description

---

Hi-Media Pseudomonas Agar For Detection of Pyocyanin, 500 g

## Additional Information

---

No.	MUI19-500G
Manufacturer (Brand)	Hi-Media
HSNumber	38210000
Transport temperature	Room temperature
autoclavable at 121°C	No
Pack	500 g
Storage	Room temperature
sterile	No

