



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



HiMedia Reinforced Clostridial Agar, granulated, 500 g

€81.50
plus VAT &
Shipping

Product Images



Description

Hi-Media Reinforced Clostridial Agar, granulated, 500 g

Additional Information

| | |
|-----------------------|------------------|
| No. | GM154-500G |
| Manufacturer (Brand) | Hi-Media |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Pack | 500 g |
| Storage | Room temperature |
| sterile | No |

