



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoFroxx Potassium Chloride - Solution 3 M, saturated with silver chloride for storing pH-

**€32.00** As low

as

**€28.25**

**plus VAT &  
Shipping**

### Product Images

## Description

---

Trust is good, control is better! To regularly check the accuracy of your pH electrode – and readjust it if necessary – we offer suitable buffer solutions with precisely defined pH values. The pH calibration solutions are available as colorless and colored versions (red, green, blue). The batch-specific certificate of analysis guarantees a maximum deviation of  $\pm 0.02$  pH units.

## Additional Information

---

No.	LC-7035.2
Manufacturer (Brand)	neoFroxx
EAN	4066306032981
HSNumber	38220000
Transport temperature	Room temperature
autoclavable at 121°C	No
Pack	500 ml
Storage	Room temperature
sterile	No

