



neoLab® neoLab 19TLG. Pliers/screwdriver set

Product Images



Description

neoLab Toolbox 19 parts - Complete tool set

Scope of delivery

- Stainless steel tweezers
- 12 computer screwdrivers: slotted 2/2.5/3 mm, PH 00/0/1, T-profile T6/T7/T8/T9/T10/T15 mm
- 6 mini pliers: Side cutter 115 mm, needle nose pliers straight 115/125 mm, needle nose pliers curved 125 mm, snipe nose pliers 115 mm, flat nose pliers 125 mm
- Black bag

Additional Information

No.	DK-0565
Manufacturer (Brand)	neoLab
EAN	4066292505995
sterile	No

