



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® box, with lid, 60 x 40 mm

**€61.39** As low as  
**€47.22**  
plus VAT &  
Shipping

## Product Images



## Description

DURAN® can, with lid, Ø 60 x 40 mm For storing substances and materials, it is ideal to use containers that can be closed. In this way, substances and materials can be protected from undesirable effects. These vessels are available in different shapes such as cans, cylinders or boxes. The choice of vessel is made according to the application and/or nature of the material. Special features and advantages: With lid for safe storage and transport Baked and therefore permanently durable white printed image Application: Storage and transport

## Additional Information

No.	SD-1243
Manufacturer (Brand)	DURAN
VGKL number	122420834
Transport temperature	Room temperature

