



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Polycarbonate Straight-Walled Wide-Mouth Containers with Closure

€21.68 As low as
€17.52
**plus VAT &
Shipping**

Product Images



Description

Nalgene Polycarbonate Straight-Walled Wide-Mouth Containers with Closure Use these clear polycarbonate straight-walled wide-mouth containers to collect and store pastes, solids, samples and lab supplies: highly impact resistant and extremely versatile. Capacity: 60 ml, 2 oz. Capacity (metric): 60 ml capacity (anglo): 2 oz. Inner Diameter (Anglo) Neck: 1.7 in. Inner Diameter (metric) Neck: 43 mm For use with (Application): Storage Contains: Cover Quantity per box: 48 Quantity per package: 4 Type: Straight Sided Jar Outside diameter (anglo): 1.9 in. Outer Diameter (metric): 48 mm Clarity: Clear Closure Size: 53 mm Closure Material: Polypropylene Closure Type: Screw Height (Anglo): 1.8 in. Height (metric): 45 mm Opening: Wide Shape: Round Unit Size: Case of 48

Additional Information

No.	CD-0388
Manufacturer (Brand)	Nalgene
VGKL number	228721162
Transport temperature	Room temperature

