



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Polycarbonate Straight-Walled Wide-Mouth Containers with Closure

€186.96
plus VAT &
Shipping

Product Images



Description

Nalgene Polycarbonate Straight-Walled Wide-Mouth Containers with Closure Use these clear polycarbonate straight-walled wide-mouth containers to collect and store pastes, solids, samples and lab supplies: highly impact resistant and extremely versatile. Capacity: 125 ml (4 oz.) capacity (metric): 125 ml capacity (anglo): 4 oz. Inner diameter (anglo) Neck: 2.5 in. Inner diameter (metric) Neck: 64 mm For use with (application): Storage Contains: Cover Quantity per carton: 24 Quantity per package: 4 Type: Straight-walled vessel Outer diameter (Anglo-American): 2.5 in. Outer diameter (metric): 64 mm Clarity: clear Closure size: 70 mm Closure material: polypropylene Closure type: screw Height (anglo): 2.9 in. Height (metric): 74 mm Opening: Wide Shape: Round Unit Size: Case of 24

Additional Information

No.	CD-0389
Manufacturer (Brand)	Nalgene
VGKL number	228721163
Transport temperature	Room temperature

