




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® straight-walled wide-mouth vessels made of PMP with white

€79.76 As low
as
€64.24
plus VAT &
Shipping

Product Images



Description

Nalgene PMP straight wall wide mouth jars with white polypropylene screw cap Protect against breakage and contamination with these museum quality PMP jars that lose only 2% formaldehyde per year. Capacity: 1000 ml, 32 oz. Capacity (metric): 1000 ml capacity (anglo): 32 oz. Inside diameter (anglo) neck: 4.4 in. Inner Diameter (metric) Neck: 112 mm Includes: Screw Closure Quantity Per Case: 16 Quantity Per Package: 4 Type: Jar Outer Diameter (Anglo): 4.4 in. Outer Diameter (metric): 112 mm Clarity: Clear Closure Size: 120 mm Closure Material: polypropylene Closure Type: screw Height (anglo): 5.9 in. Height (metric): 151 mm Opening: Wide Shape: Round Unit Size: Case of 16

Additional Information

No.	CD-0397
Manufacturer (Brand)	Nalgene
VGKL number	228721176
Transport temperature	Room temperature

