



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene® straight-walled wide-mouth vessels made of PMP with white

**€42.20** As low  
as  
**€34.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Nalgene PMP straight wall wide mouth jars with white polypropylene screw cap Protect against breakage and contamination with these museum quality PMP jars that lose only 2% formaldehyde per year. 250 ml (8 oz.) capacity (metric): 250 ml Capacity (anglo): 8 oz. Inner Diameter (Anglo) Neck: 2.5 in. Inner Diameter (metric) Neck: 64 mm Includes: Screw Closure Quantity Per Case: 24 Quantity Per Package: 4 Type: Jar Outer Diameter (Anglo): 2.9 in. Outer Diameter (metric): 75 mm Clarity: Clear Closure Size: 70 mm Closure Material: Polypropylene Closure Type: Screw Height (Anglo): 4.7 in. Height (metric): 119 mm Opening: Wide Shape: Round Unit Size: Case of 24

## Additional Information

---

No.	CD-0395
Manufacturer (Brand)	Nalgene
VGKL number	228721174
Transport temperature	Room temperature

