




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Heavy Duty Polypropylene Vacuum Balloon Bottle with

€639.00
plus VAT &
Shipping

Product Images

Description

Heavy Duty Nalgene Polypropylene Vacuum Balloon Bottle with Handle These heavy duty vacuum balloon bottles are suitable for the most extreme operating conditions, can be used as vacuum traps and hold a full vacuum for up to 8 hours. Capacity (Anglo-American): 5 gal. Closure Material: polypropylene Closure Size: 83 B Description: heavy duty vacuum balloon bottles Outer Diameter (metric): 282 mm Capacity (metric): 20 l Type: balloon bottle with handle For use with (application): Ideal for the most extreme operating conditions Contains: Handles Quantity per carton: 4 Quantity per package: 1 Volume to rim (anglo): 5.5 gal. Volume to rim (metric): 22.5 l Autoclavable: Yes Unit Size: Case of 4

Additional Information

No.	CD-0533
Manufacturer (Brand)	Nalgene
VGKL number	228722260
Transport temperature	Room temperature

