



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® HDPE, Amber Carboy

€567.00
plus VAT &
Shipping

Product Images

Description

Nalgene HDPE, Amber Carboy Can be used with: To protect your light-sensitive chemicals, reagents, buffers and standards during mixing and storage; Closure material: polypropylene; Color: Amber Capacity (Anglo-American): 2.5 gal. Certification/Compliance: meets light transmission test requirements (US Pharmacopeia latest edition) for storage of light sensitive substances. Cap Material: Polypropylene Cap Size: 83B Description: Brown HDPE Balloon Bottle Capacity (metric): 10 l Type: Carboy For use with (application): Protect photosensitive chemicals, reagents, buffers and standards when mixing and storing. Includes: Wide shoulder handles Quantity per carton: 6 Quantity per package: 6 Volume to rim (anglo): 3.1 gal. Volume to rim (metric): 12.5 l Unit Size: Case of 6

Additional Information

No.	CD-0554
Manufacturer (Brand)	Nalgene
VGKL number	228722567
Transport temperature	Room temperature

