



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® round LDPE balloon bottles with tap

€209.00
plus VAT &
Shipping

Product Images



Description

Nalgene round LDPE balloon bottles with stopcock Use the Thermo Scientific Nalgene round LDPE balloon bottles with stopcock to collect and dispense distilled water, acids and reagents. Capacity (Anglo-American): 13 gal. Cap material: polypropylene, thermoplastic elastomer seal Cap size: 83 B Description: round balloon bottles with stopcock Inner diameter (metric) Neck: 64 mm Type: balloon bottle with stopcock For use with (application): For collecting and dispensing distilled water, acids and reagents Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with cap: 678 mm Includes: Integrated shoulder handles, molded-in stopcock Autoclavable: Non-autoclavable Color: White (opaque) Unit Size: Each

Additional Information

No.	CD-0568
Manufacturer (Brand)	Nalgene
VGKL number	228723185
Transport temperature	Room temperature

