



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® LDPE, Wide-Mouth Carboy with Handle

€315.00
plus VAT &
Shipping

Product Images

Description

Nalgene LDPE, Wide-Mouth Carboy with Handle The Thermo Scientific Nalgene wide-mouth balloon bottles made of LDPE with handles are suitable for storing and transporting solids and powders. Capacity (Anglo-American): 2.5 gal. Closure Material: Polypropylene Closure Size: 100 - 415 mm Capacity (metric): 10 l Type: Carboy with Handle For use with (application): Store and transport solids, powders and liquids Number per case: 6 Number per package: 1 Volume to rim (Anglo-American): 3 gal. Volume to rim (metric): 12 l Autoclavable: Not Autoclavable Unit Size: Case of 6

Additional Information

No.	CD-0536
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

