



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene® Clearboy, round, polycarbonate with tap, 10 l

**€969.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Nalgene Clearboy, round, polycarbonate with tap, 10 l Thermo Scientific Nalgene Clearboys with tap are durable, lighter and safer than glass. Capacity (Anglo-American): 5 gallon. Closure material: polypropylene, thermoplastic elastomer seal Closure size: 83B Description: transparent balloon bottle with stopcock made of PC, 20 l inner diameter (metric) Neck: 66 mm Type: balloon bottle with stopcock For use with (application): Suitable for mixing, storing and freezing aqueous reagents and biologicals, especially useful for visual inspection of contents for quality. Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with cap: 536 mm Contains: Tap Number per case: 4 Number per package: 1 Autoclavable: Yes Unit Size: Case of 4

## Additional Information

---

No.	CD-0563
Manufacturer (Brand)	Nalgene
VGKL number	228723172
Transport temperature	Room temperature

