



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Polypropylene, Carboy with Sanitary Flange

€548.00
plus VAT &
Shipping

Product Images

Description

Nalgene Polypropylene, Carboy with Sanitary Flange Use these balloon bottles for dispensing or, with an appropriate connector, as reagent containers to supply systems such as bioreactors or chromatography columns. Capacity (Anglo-American): 2.6 gal. Certification/Compliance: USP Class VI, non-cytotoxic Closure Material: Polypropylene Closure Size: 83 B Description: Autoclavable Carboy Sanitary Flange Capacity (metric): 10 L Type: Carboy with Sanitary Flange For use with (application): Supply reservoirs to systems such as fermentors or chromatography columns. Graduations: Graduated to contain in liters and gallons Contains: Molded-in handles, 1-1/2 in. secure, non-threaded sanitary dispensing port Volume to rim (Anglo-American): 3 gal. Volume to rim (metric): 12 L Autoclavable Unit Size: Each

Additional Information

No.	CD-0674
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

