



Nalgene® Polycarbonate Biotainer bottles and

Product Images



Description

Nalgene Polycarbonate Biotainer Bottles and Balloon Vials Reduce the risk of carryover contamination and save the cost of cleaning, sterilization, and related validations. Use sterile, ready-to-use Thermo Scientific Nalgene Biotainer Bottles and Balloon Bottles, designed specifically for pharmaceutical and biological product manufacturers. Excellent for freezing. Capacity (metric): 20 mL Capacity (anglo): 0.7 oz. Inner Diameter (metric) Neck: 20mm Quantity Per Case: 500 Quantity Per Package: 100 Type: bottle and balloon bottle Lining Material: silicone Certification/Compliance: USP Class VI, non-cytotoxic, USP It;661 gt;, non-pyrogenic, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 20 mm Closure Material: Polypropylene For use with (devices): Freezing Graduated: Yes Height (Anglo-American): 2.7 in. Opening: Narrow Shape: Round Unit Size: Case of 500

Additional Information

No.	CD-1204
Manufacturer (Brand)	Nalgene

