




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene® autoclavable polypropylene balloon bottle with tap, 10 l

**€662.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Nalgene autoclavable polypropylene balloon bottle with stopcock, 10 l Thermo Scientific Nalgene autoclavable clarification bottles with stopcock have excellent chemical resistance and are suitable for distilled water as well as nutrient media. Capacity (Anglo-American): 5.25 gal . Closure material: polypropylene, thermoplastic elastomer seal Closure size: 83B Description: autoclavable balloon bottle with stopcock Inner diameter (metric) Neck: 64 mm Outer diameter (metric): 286 mm Capacity (metric): 20 l Type: balloon bottle with stopcock For use with (application): Solutions and media for long-term storage and dispensing Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with cap: 528 mm Includes: Integral shoulder handles, molded-in stopcock Number per case: 4 Number per package: 1 Autoclavable: Yes Unit Size: Case of 4

## Additional Information

---

No.	CD-0570
Manufacturer (Brand)	Nalgene
VGKL number	228723192
Transport temperature	Room temperature

