



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Clearboy, round, polycarbonate with tap, 10 l

€747.00
plus VAT &
Shipping

Product Images



Description

Nalgene Clearboy, round, polycarbonate with tap, 10 l Thermo Scientific Nalgene Clearboys with tap are durable, lighter and safer than glass. Capacity (Anglo-American): 2.5 gal. Closure Material: polypropylene, thermoplastic elastomer seal Closure Size: 83B Description: transparent balloon bottle with stopcock made of PC, 10 l inner diameter (metric) Neck: 66 mm Type: balloon bottle with stopcock For use with (application): Suitable for mixing, storing and freezing aqueous reagents and biologicals, especially useful for visual inspection of contents for quality. Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with cap: 394 mm Contains: Tap Number per case: 4 Number per pack: 1 Autoclavable: Yes Unit Size: Case of 4

Additional Information

No.	CD-0562
Manufacturer (Brand)	Nalgene
VGKL number	228723171
Transport temperature	Room temperature

