




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Polypropylene, Wide-Mouth Carboy with Handle

€539.00
plus VAT &
Shipping

Product Images

Description

Nalgene Polypropylene, Wide-Mouth Carboy with Handle The convenient wide mouth of this balloon bottle makes cleaning and dispensing easy. The wide-mouthed handles make it easy to carry and pour, even with gloves. Capacity (Anglo-American): 5.3 gal. Closure Material: polypropylene Closure Size: 100-415 mm Description: autoclavable wide-neck balloon bottle with handles Outer Diameter (metric): 282 mm Capacity (metric): 20 l Type: Carboy with Handle For use with (Application): Store and transport solids, powders and liquids Graduations: Graduation in liters Height (metric) with closure: 483 mm Contains: Handles Quantity per case: 4 Quantity per package: 1 Volume to rim (Anglo): 5.5 gal. Volume to rim (metric): 22.5 l Autoclavable: Yes Mouth Style: Wide Unit Size: Case of 4

Additional Information

No.	CD-0540
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

