



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene® Polypropylene, Wide-Mouth Carboy with Handle

**€94.60** As low  
as  
**€76.50**  
plus VAT &  
Shipping

### Product Images

## Description

---

Nalgene Polypropylene, Wide-Mouth Carboy with Handle The convenient wide mouth of this balloon bottle makes cleaning and dispensing easy. The wide-mouthed handles make it easy to carry and pour, even with gloves. Capacity (Anglo): 2.5 gal. Closure Material: polypropylene Closure Size: 100 - 415 mm Description: autoclavable wide-neck balloon bottle with handles Outer Diameter (metric): 250 mm Capacity (metric): 10 l Type: Carboy with Handle For use with (Application): Store and transport solids, powders and liquids Graduations: Graduation in liters Height (metric) with closure: 343 mm Contains: Handles Quantity per case: 6 Quantity per package: 1 Volume to rim (Anglo-American): 3 gal. Volume to rim (metric): 12 l Autoclavable: Yes Mouth Style: Wide Unit Size: Case of 6

## Additional Information

---

No.	CD-0539
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

