



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® Certified Clean Polycarbonate Biotainer Carboys

€712.02
plus VAT &
Shipping

Product Images

Description

Nalgene Certified Clean Polycarbonate Biotainer Carboys Specifically designed for pharmaceutical and biological drug manufacturing to reduce the risk of carryover contamination. Durable translucent PC construction ensures sample integrity, even when freezing to -100 C. Capacity (metric): 5 L Capacity (Anglo): 169.1 oz. Color: Blue-Green Inner Diameter (metric) Neck: 48 mm Number Per Case: 6 Number Per Package: 1 Type: Biotainer Bottles and Balloon Bottle Closure Color: White Lining Material: Silicone Certification/Compliance: USP Class VI, non-cytotoxic, USP It;661gt;, non-pyrogenic, FDA 21CFR 177.1315 Closure Included: Yes Autoclavable: Not autoclavable Closure Size: 48 mm Closure Material: Polypropylene Diameter (metric): 166 mm Graduated: Yes Height (Anglo-American): 11.8 in. Height (metric): 299 mm Opening: One Narrow Shape: Square Unit Size: Case of,6

Additional Information

No.	CD-1396
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

