



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle® storage bottle o. threaded spout, HDPE, 25 l, m.V.

€95.75 As low
as
€82.60
plus VAT &
Shipping



Description

Extra thick-walled, sturdy storage bottle with sturdy carrying handle for safe transport and storage. Printed with blue liter scale. The aeration/deaeration valve (accessory) allows air to enter the container to avoid a vacuum in the bottle. The canister can also be used very well as a clarification canister. A drain cock or other accessories with 3/4" female thread/ union nut can be connected via a threaded connection. HDPE bottle opening diameter 50 mm

Additional Information

No.	KL-2011
Manufacturer (Brand)	Bürkle
VGKL number	220530203
Transport temperature	Room temperature

