



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® LDPE, round balloon bottle

€337.98
plus VAT &
Shipping

Product Images

Description

Nalgene LDPE, round balloon bottle Thermo Scientific Nalgene LDPE balloon bottles with handles are the ideal solution for storing and transporting reagents. Capacity (Anglo-American): 2.5 gal. Cap Material: Polypropylene Cap Size: 83B Description: Round LDPE Balloon Bottle with Handle Inner Diameter (metric) Neck: 64 mm Outer Diameter (metric): 250 mm Capacity (metric): 10 l Type: Carboy with Handle For use with (application): Store and transport reagents Height (metric) with closure: 389 mm Contains: Handles Number per box: 6 Number per package: 1 Autoclavable: Not Autoclavable Unit Size: Case of 6

Additional Information

No.	CD-0512
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

