



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® LDPE, round balloon bottle

€415.00
plus VAT &
Shipping

Product Images

Description

Nalgene LDPE, round balloon bottle Thermo Scientific Nalgene LDPE balloon bottles with handle are the ideal solution for storing and transporting reagents. Capacity (Anglo-American): 6.5 gal . Cap Material: polypropylene Cap Size: 83B Description: round LDPE balloon bottle with handle Inner Diameter (metric) Neck: 64 mm Outer Diameter (metric): 287 mm Capacity (metric): 25 l Type: balloon bottle with handle For use with (application): Storage and transport of reagents Height (metric) with closure: 594 mm Includes: Handles Quantity per box: 4 Quantity per package: 1 Autoclavable: Not autoclavable Unit Size: Case of 4

Additional Information

No.	CD-0515
Manufacturer (Brand)	Nalgene
Transport temperature	Room temperature

