



Nalgene® Autoclavable Measuring Bucket

Product Images



Description

Nalgene Autoclavable Graduated Measuring Buckets The durable and convenient Thermo Scientific™ Nalgene™ Autoclavable Graduated Laboratory Buckets allow you to store a variety of substances, including chemical and corrosive liquids. Description: 20 x 20cm ; Instrument only Autoclavable Handle: Yes Capacity (anglo): 8 qt. Graduations: Molded Item Description: Bucket Graduated PP

Additional Information

No.	CD-1803
Manufacturer (Brand)	Nalgene

