



Nalgene® Autoclavable Measuring Bucket

Product Images



Description

Nalgene Autoclavable Measuring Buckets Store a variety of substances, including chemical and corrosive liquids, in the durable and convenient Thermo Scientific™ Nalgene™ Autoclavable Graduated Laboratory Buckets. Description: graduated bucket, 14 qt. Autoclavable Handle: Yes Capacity (Anglo-American): 14 qt. Graduations: Molded Item Description: Bucket Graduated PP 13.2 L

Quantity Per Package: 1 Quantity Per Case: 6 Capacity: 14 qt. (13.2 liters) Unit Size: Case of 6

Additional Information

No.	CD-1805
Manufacturer (Brand)	Nalgene

