



Nalgene® LDPE airtight buckets

Product Images

Description

Nalgene Airtight LDPE Pails Thermo Scientific Nalgene pails are made of LDPE and meet the US OSHA standard for hematogenous pathogens with their airtight lids. Autoclavable: Autoclavable Handle: Yes Capacity (anglo): 8 qt. Capacity (metric): 7.6L Graduations: Molded Height (Anglo-American) With Closure: In 10.5 . Height (metric) with closure: 26.7cm Item Description: Pail Air-Tight LDPE 7.6 L Lid: Yes Certification/Compliance: Meet OSHA 29CFR 1910.1030 standard for use as protection against bloodborne pathogens. Quantity per Package: 1 Quantity per Case: 12 Capacity: 8 qt., 7.6L Unit Size: Case of 12

Additional Information

No.	CD-1806
Manufacturer (Brand)	Nalgene

