



Nalgene® HDPE rectangular balloon bottle

Product Images

Description

Nalgene HDPE, rectangular balloon bottle The wide opening of these rectangular balloon bottles made of HDPE simplifies filling, transfer and cleaning. They impress with their sturdy, space-saving design. Capacity (Anglo-American): 2.4 gal. Closure Material: polypropylene Closure Size: 100 - 415 mm Capacity (metric): 9 L Dimensions (L x W x H): 8.5 x 5.75 x 13.5 in. (216 x 146 x 343 mm) Type: Rectangular Balloon Bottle For use with (Application): Convenient for dry materials such as solids and powders Height (Anglo-American): 13.5 in. Height (metric): 343 mm Length (anglo): 8.5 in. Length (metric): 216 mm Quantity per carton: 6 Quantity per package: 1 Width (anglo): 5.75 in. Width (metric): 146 mm Autoclavable: Not autoclavable Unit Size: Case of 6

Additional Information

No.	CD-0521
Manufacturer (Brand)	Nalgene

