



# Nalgene® HDPE rectangular balloon bottle

## Product Images

## Description

Nalgene HDPE, rectangular balloon bottle The wide opening of these rectangular balloon bottles made of HDPE simplifies filling, transfer and cleaning. They feature a rugged, space-saving design. Capacity (Anglo-American): 5 gal. Closure material: polypropylene Closure size: 100 - 415 mm Capacity (metric): 20 l Dimensions (L x W x H): 12.5 x 8.63 x 15 in. (318 x 219 x 381 mm) Type: Rectangular Carboy For use with (application): Convenient for dry materials such as solids and powders Height (Anglo-American): 15 in. Height (metric): 381 mm Length (anglo): 12.5 in. Length (metric): 318 mm Number per carton: 4 Number per package: 1 Width (anglo): 8.63 in. Width (metric): 219 mm Autoclavable: Not Autoclavable Unit Size: Case of 4

## Additional Information

|                      |         |
|----------------------|---------|
| No.                  | CD-0522 |
| Manufacturer (Brand) | Nalgene |

