



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Linker® wire baskets made of welded mesh w 8 mm, d 1 mm, for e.g. test tubes Nominal size: 42 x 20 x 20 cm

€102.00 As low
as
€88.00
plus VAT &
Shipping

Product Images



Description

Wire baskets, made of welded mesh w 8 mm, d 1 mm, for e.g. test tubes etc., with white plastic cover, nominal size: 42 x 20 x 20 cm

Additional Information

| | |
|-----------------------|------------------|
| No. | LR-0011 |
| Manufacturer (Brand) | Linker |
| Transport temperature | Room temperature |

