




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Evaporating dish, with spout, 600ml

€282.50
plus VAT &
Shipping

Product Images



Description

DURAN® Evaporating dish, 600 ml Converting substances and liquids into other states of aggregation and thus initiating further reactions is very important for daily work in the laboratory. A wide variety of vessels, apparatus and equipment are used for these activities. Special features and advantages: Precise spout for clean work Easy-to-read scale with line graduation for measuring the contents Baked-on and therefore permanently durable white printed image Application: Evaporation of liquids, preparation of concentrated solutions, prevention of solvent condensation on glass walls, sizes 25, 30, 35 and 45 ml without title block.

Additional Information

| | |
|-----------------------|------------------|
| No. | SD-0447 |
| Manufacturer (Brand) | DURAN |
| VGKL number | 112130154 |
| Transport temperature | Room temperature |

