



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Crystallizing dish, without spout, 3500 ml

€129.70 As low
as
€99.77
plus VAT &
Shipping

Product Images



Description

DURAN® Crystallizing dish, without spout, 3500 ml Converting substances and liquids into other states of aggregation and thus initiating further reactions is very important for daily work in the laboratory. A wide variety of vessels, apparatus and equipment are used for these activities. Special features and advantages: Flat and wide shape enable fast and uniform evaporation of solvents Large base area Precise spout for clean work Burnt-in and therefore permanently durable white print Application: crystallization , preparation of concentrated solutions, Available with spout DIN 12338 and without spout DIN 12337.

Additional Information

No.	SD-0469
Manufacturer (Brand)	DURAN
VGKL number	112131363
Transport temperature	Room temperature

