



Nalgene® autoclavable polypropylene trays

Product Images

Description

Nalgene™ Autoclavable Polypropylene Trays Use the durable Thermo Scientific™ Nalgene™ autoclavable polypropylene trays to repeatedly autoclave instruments and laboratory equipment. Dimensions (L x W x H): 10.25 x 6.25 x 4.125 in. (260 x 159 x 105 mm) Description: Autoclavable polypropylene tray Type: tray Length (Anglo-American): 10.25 in. Quantity per package: 1 Height (Anglo-American): 4.125 in. Shape: Rectangular Width (Anglo): 6.25 in. Autoclavable: Yes Unit Size: Case of 6

Additional Information

No.	CD-1794
Manufacturer (Brand)	Nalgene

